



MedGyn

Innovative Products for  
Ob/Gyn Professionals

## MedGyn (Office Use) Endoscopic (Hysteroscopic) Cannula

*MedGyn's Endoscopic (Hysteroscopic) Cannula is a single-use cannula for uterine observation pre- and post-surgery.*



### Features & Diagnostic Benefits

- Sterile, single-use, highly-durable lightweight, water resistant 6mm, 7mm, 8mm, and 9mm cannulas
- High resolution eagle eye optical system
- Provides viewing of entire uterine cavity
- Anatomical curve designed cannula with dual function: Fluid infusion and suction
- Cannula has connections for light source and 8" monitor
- Ideal diagnostic device for: Cases of abnormal uterine bleeding, infertility, uterine polyps and fibroids, suspected uterine or cervical cancer, repeated miscarriages and uterine malformations
- Workstation includes cannula, monitor and foldable stand. Up-right rolling stand is optional



Monitor with foldable Stand



Endoscopic Cannula (Disposable)



# Therapeutic For

- Removal of endometrial polyps
- Removal of small fibroids
- Removal of adhesions (Asherman's Syndrome)
- IUD setup and removal
- Dilation and Curettage

## Specifications

Illumination > 500 LX

Field of View > 60 Degree

Resolution > 201 p/mm

Working Distance 170 mm



Monitor with up-right rolling Stand

- 018145** MedGyn Endoscopic (Hysteroscopic) Cannula
- 018146** MedGyn Endoscopic (Hysteroscopic) Cannula Monitor with up-right rollind Stand
- 018147** MedGyn Endoscopic (Hysteroscopic) Cannula Monitor with foldable Stand

## TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

### MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

**t:** +1 630.627.4105 toll-free: 800.451.9667

**f:** +1 630.627.0127

**e:** info@medgyn.com

**w:** medgyn.com

