


Cell Sweep®

The MedGyn Cell Sweep® is an all-in-one pap smear device that eliminates the need for a cytobrush and spatula. Cell Sweep® completely embraces and covers the T-zone/SCJ during cell sampling.



Features & Benefits

- Cell Sweep® allows for optimal cell harvest and cell transferring to achieve an effective pap smear
- Dual sampling of Endo and Ectocervi
- Thin layer preparation of slide using liquid base transfer medium
- Separation of endocervical and ectocervical cells in a one slide technique
- Specific design ensures achieving a consistent monolayered spread
- Minimizes false negatives as it completely covers the T-zone/SCJ
- Sexually transmitted disease sampling
- Minimizes red blood cell and blood artifact



Simplicity

Cell Sweep® is a single device with the capacity for simultaneous endocervical and ectocervical sampling. This is normally achieved using a brush and a spatula separately. It is soft, malleable, and adaptable to individual cervical variations. It samples the entire T-zone while remaining in close contact with the cervix.

The design and specific structure of Cell Sweep® enables the cervix to be sampled twice, simultaneously, in a 360° motion. While the bristles are harvesting endocervical cells, both the spatula section and four lower, longer bristles are harvesting ectocervical cells.



- 022363 MedGyn Cell Sweep® Sterile (100/pk)
- 022361 MedGyn Cell Sweep® Non Sterile (100/pk)
- 022363K MedGyn Cell Sweep® Sterile with Kit (Slide and Slide Holder) (100/pk)
- 022361K MedGyn Cell Sweep® Non Sterile with Kit (Slide and Slide Holder) (100/pk)
- 022362 MedGyn Cell Sweep® Slide Holder (100/pk)
- 022361K MedGyn Cell Sweep Non Sterile with Kit (Slide and Slide Holder)
- 022363K MedGyn Cell Sweep Sterile with Kit (Slide and Slide Holder)

TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

t: +1 630.627.4105 toll-free: 800.451.9667

f: +1 630.627.0127

e: info@medgyn.com

w: medgyn.com

