

MedGyn Disposable Electrodes

The MedGyn's Disposable Electrodes are made with high quality tungsten wire for superior accuracy. MedGyn offers a full range of tips and sizes to suit most electrosurgical procedures. Because it retains its mechanical properties at higher temperature than stainless steel, tungsten maintains better shape and wire integrity throughout the excision procedure.

Features & Benefits

- Complete selection includes loops, balls, and needle tips specifically designed for electrosurgical excision procedures
- Individually-wrapped and sterile; eliminates cross-contamination
- Utilizes standard 2.3mm stainless steel shaft



Ball Electrodes

028041 Ball Electrode 3mm Ball(/pk)

028042 Ball Electrode 5mm Ball(5/pk)


Square Electrode

028075 Square Loop 5mm x 5mm (5/pk)

028072 Square Loop 10mm x 4mm (5/pk)

028073 Square Loop 10mm x 8mm (5/pk)

028074 Square Loop 10mm x 10mm (5/pk)



Blade Electrodes

028081 Blade Electrode (5/pk)

028053 Needle 5mm (5/pk)

Loop Electrodes

028067 Loop 5mm x 5mm (5/pk)

028069 Loop 10mm x 5mm (5/pk)

028070 Loop 10mm x 7mm (5/pk)

028071 Loop 10mm x 8mm (5/pk)

028055 Loop 10mm x 10mm (5/pk)

028056 Loop 10mm x 13mm (5/pk)

028057 Loop 15mm x 5mm (5/pk)

028087 Loop 15mm x 7mm (5/pk)

028058 Loop 15mm x 8mm (5/pk)

028068 Loop 15mm x 10mm (5/pk)

028076 Loop 15mm x 13mm (5/pk)

028059 Loop 15mm x 15mm (5/pk)

028060 Loop 20mm x 8mm (5/pk)

028063 Loop 20mm x 10mm (5/pk)


028061 Loop 20mm x 12mm (5/pk)

028062 Loop 20mm x 13mm (5/pk)

028064 Loop 20mm x 15mm (5/pk)

028065 Loop 20mm x 20mm (5/pk)

028066 Loop 25mm x 10mm (5/pk)



TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

t: +1 630.627.4105 **toll-free:** 800.451.9667

f: +1 630.627.0127

e: info@medgyn.com

w: medgyn.com

