



MedGyn

Innovative Products for
Ob/Gyn Professionals

Pipette™, Pipette II™ & IV™

MedGyn's Pipette™, Pipette II™ & Pipette IV™ are unique, sterile, disposable cannula for sampling the uterine endometrium.

Simplicity

The MedGyn Pipette™, Pipette II™, and Pipette IV™ is utilized to obtain clearly differentiated endometrial tissue without anesthesia, in an outpatient setting or as an office procedure. This device is designed to improve the ease of tissue collection and to reduce the physical discomfort to the patient, which is inherent to a biopsy.

It provides sufficient histologic specimen and is procedurally cost-effective for routine use. A cervical block, dilator or tenaculum is not required for using this device.

The MedGyn Pipette™, Pipette II™, Pipette IV™ is recommended for both pre- and post-menopausal women for diagnostic purposes.



Pipette



Pipette II



Pipette IV





MedGyn

Innovative Products for
Ob/Gyn Professionals

Features & Benefits

- Convenient for routine gynecological use.
- Facilitates endometrial aspiration for specimen removal.
- Helps avoid endometrial sample 'contamination' with endocervical cells by collecting sample prior to curette removal.
- Semi-rigid tapered 3.1 mm cannula enables easier insertion and minimizes the possibility of perforation.
- Yield specimens highly suitable for both cytologic and histologic diagnosis.
- Flexible tube bent conveniently with sharp slot on side assures better sampling from all areas of the endometrial cavity, even from the most difficult to reach areas.

MedGyn Pipette's

022721 MedGyn Pipette 3.1mm O.D - One hole (25/bx)

022734 MedGyn Pipette II 3.1 mm O.D - Two holes (25/bx)

022730 MedGyn Pipette IV 3.1 mm O.D - Four holes (25/bx)

Trusted Products for Premier Women's Care

MedGyn Products Inc.

100 W. Industrial Rd. Addison, IL 60101 USA

t: 630.627.4105 toll-free: 800.451.9667

f: 630.627. 0127

e: info@medgyn.com

w: medgyn.com

