

MedGyn Laminaria “Biodegradable”

MedGyn's Laminaria is an all-natural, disposable and sterile cervical dilator made from a sea-grown plant containing no synthetic materials. Laminaria provides gentle, gradual dilation.

Features & Benefits

MedGyn's disposable and sterile Laminaria is an all-natural cervical dilator made from a sea-grown plant containing no synthetic materials. Safely used for over a century, Laminaria provides gentle, gradual dilation. Once inserted into the cervical canal, Laminaria absorbs liquid, swelling several times its initial dried diameter while gently dilating the cervix.



Medgyn's McMahon
Laminaria Inserter

MedGyn's Laminaria
Insertor



Indication for use:

Laminaria “Biodegradable” should be considered whenever cervical dilation and/or softening of the cervix is desired. Some examples included: cervical stenosis, dysmenorrhea, placement and removal of intrauterine devices, induction of labor, radium placement, curettage and drainage of uterine cavity, endometrial biopsy, and suction cannula aspiration. is slightly unclear and doesn't encompass the full spectrum of uses. Also, a lot of the physicians know it's used for D&C, so we should increase the importance of the other uses. Here is a revised version:

- Gradual Induction of Labor
- Diminution of risk of uterine perforation during hysteroscopy
- Alleviation of cervical stenosis
- Ease of insertion and removal of intrauterine devices
- Facilitation of insertion of radium
- Minimizing pain of cervical dilation prior to endometrial biopsy or uterine curettage/suction aspiration

With Button (prevents over-insertion)

021022 Laminaria #2 Extra Small w/button

021023 Laminaria #3 Small w/button

021024 Laminaria #4 Medium w/button

021025 Laminaria #5 Large w/button

021026 Laminaria #6 Extra Large w/button

Laminaria Inserter

031201 MedGyn's Laminaria Inserter

031202 MedGyn's McMahon Laminaria Inserter



MedGyn's Laminaria
Inserter



MedGyn's McMahon
Laminaria Inserter

Standard, approx. length 60mm

021002 Laminaria #2 Small Thin

021003 Laminaria #3 Small

021004 Laminaria #4 Medium

021005 Laminaria #5 Large/Thick

021006 Laminaria #6 Extra Large Thick

021008 Laminaria #8 Jumb

021010 Laminaria #10 Extra Jumbo

Extra Long, approx. 70 mm

021014 Laminaria #4, Long(70mm)

021015 Laminaria #5, Long(70mm)

021016 Laminaria #6, Long(70mm)

021018 Laminaria #8, Long(70mm)

TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

t: +1 630.627.4105 toll-free: 800.451.9667

f: +1 630.627.0127

e: info@medgyn.com

w: medgyn.com

